

# Have you ever had an ICU patient...

- ▶ become progressively more swollen and edematous after fluid resuscitation?
- ▶ develop renal failure and need dialysis?
- ▶ develop multiple organ failure and die?

## What was their Intra-Abdominal Pressure?

AbViser® AutoValve® IAP Monitoring Device is a closed IAP monitoring system with proven accuracy and reproducibility.<sup>1</sup>

- **Simplifies ICU Nursing Tasks** –  
Easy set-up and use
- **Standardized Data** –  
Proven accuracy and reproducibility<sup>1</sup>
- **Automatic Valve with Time Release** –  
Eliminates clamping and turning of stopcocks
- **White Stripe Tubing**  
(designates IAP) –  
Reduces misconnection errors
- **Closed System w/Syringe Sleeve** –  
Helps reduce risk of CAUTIs<sup>2</sup>
- **OR Compatible** –  
Data assists surgeon in closing laparotomy
- **Early IAH Recognition** –  
Consistent monitoring improves outcomes<sup>6,7</sup>

\*Catheter associated urinary tract infections




**Sepsis/SIRS + Volume Resuscitation → Intra-Abdominal Hypertension<sup>3-4</sup>**



**Monitoring IAP + Intervention → Reduced Organ Failure<sup>5-10</sup>**



**AbViser®  
AutoValve®**  
Intra-Abdominal Pressure Monitoring Device



## Ordering Information:

Reorder Number	Quantity / Case	Description
ABV300	10	AbViser® AutoValve® IAP Monitoring Device, Patient Mount
ABV301	10	AbViser® AutoValve® IAP Monitoring Device, Pole Mount
ABV320	10	AbViser® AutoValve® IAP Monitoring Device with Edwards Lifesciences TruWave® Transducer, Patient Mount
ABV321	10	AbViser® AutoValve® IAP Monitoring Device with Edwards Lifesciences TruWave® Transducer, Pole Mount
ABV330	10	AbViser® AutoValve® IAP Monitoring Device with ICU Medical Transpac® Transducer, Patient Mount
ABV331	10	AbViser® AutoValve® IAP Monitoring Device with ICU Medical Transpac® Transducer, Pole Mount
ABV341	10	AbViser® AutoValve® IAP Monitoring Device (Compatible with Argon CDXIII) with Transducer, Pole Mount
ABV350	5	AbViser® Neonate Adapter

[ Does not contain natural rubber latex components ] [ Individually Packaged Sterile ]

## References:

- Kimball, E.J., et al., Reproducibility of bladder pressure measurements in critically ill patients. *Intensive Care Med*, 2007.
- Huey WY, Newton DW, Augustine SC, Verjaska BD, Mitrano FP. Microbial contamination potential of sterile disposable plastic syringes. *American Journal of Hospital Pharmacy* 1985; 42:102-5
- Cheatham, M.L., et al., Results from the International Conference of Experts on Intra-abdominal Hypertension and Abdominal Compartment Syndrome. II. Recommendations. *Intensive Care Med*, 2007.
- Malbrain, M.L.N.G., et al., Incidence and prognosis of intraabdominal hypertension in a mixed population of critically ill patients: a multiple-center epidemiological study. *Crit Care Med*, 2005. 33(2): p. 315-22.
- Joseph, D.K., et al., Decompressive laparotomy to treat intractable intracranial hypertension after traumatic brain injury. *J Trauma*, 2004. 57(4): p. 687-95.
- Sun, et al., Indwelling catheter and conservative measures in the treatment of abdominal compartment syndrome in fulminant acute pancreatitis. *World J Gastroenterol*, 2006. 12(31): p. 5068-70.
- Oda, S., et al., Management of Intra-abdominal Hypertension in Patients With Severe Acute Pancreatitis With Continuous Hemodiafiltration Using a Polymethyl Methacrylate Membrane Hemofilter. *Ther Apher Dial*, 2005. 9(4): p. 355-61.
- Ivatury, R.R., et al., Intra-abdominal hypertension after life-threatening penetrating abdominal trauma: prophylaxis, incidence, and clinical relevance to gastric mucosal pH and abdominal compartment syndrome. *J Trauma*, 1998. 44(6): p. 1016-21.
- Sugrue, M., et al., Intra-abdominal hypertension is an independent cause of postoperative renal impairment. *Arch Surg*, 1999. 134(10): p. 1082-5.
- Cheatham ML, Safcsak K. Is the evolving management of intra-abdominal hypertension and abdominal compartment syndrome improving survival? *Crit Care Med*, 2010. 38(2):402-407.

® indicates trademarks of ConvaTec Inc. unless otherwise noted. TruWave is a trademark of Edwards Lifesciences Corporation and Transpac is a trademark of ICU Medical, Inc.  
© 2012 ConvaTec Inc. AP-12579-US

To learn more, call  
**1-800-422-8811**

Mon-Fri, 8:30 AM – 7:00 PM, EST  
Email: [CIC@ConvaTec.com](mailto:CIC@ConvaTec.com)  
[www.convatec.com](http://www.convatec.com)

